

Press Release

Release date: May 2013

Reach your goals easily with TSC Swab Sampling Kits

Enviroscreen Swab Sampling Kits contain a 15cm (6 inch) flexible swab, ideal for sampling hard to reach areas on equipment and surfaces.

Technical Service Consultants (TSC) Swab Sampling Kits (aka Swab Rinse Kits) are all in one sampling kits, designed to assist food operators comply with EU regulation no. 2073/2005 criteria for food safety and process hygiene, where routine microbiological testing and general hygiene monitoring of equipment and surfaces can be used as an integral part of a food safety system including HACCP, to validate and verify food safety procedures.

The peel pouch kit includes a 15cm (6inch) swab shaft ideal for reaching into awkward areas with a choice of polyester or viscose tips and a sealed transport tube containing neutralising buffer, formulated as described in ISO 18593:2004, ideal for sampling equipment or surfaces where traces of disinfectants are possible.

The 15cm swab is perfect for quick and effortless sampling of utensils, slicers, dicers, hoppers, mixers, conveyors, floors and other exposed surfaces and equipment for the recovery of TVC, gram positives (*Listeria*, *S.aureus*) and Enterobacteriaceae (*Salmonellae* and *E.coli*) and other organisms.

Benefits include:

- Simple to use
- Easy reach swab
- Blue shafts available
- Tamper evident seals
- Fits swab racks easily
- Polyester and viscose tip options
- Choice of buffers and precisely filled



TSC is ISO 9001:2000 compliant and all Enviroscreen products are sterilised by gamma irradiation and certificates of conformity can be provided for each batch. For further details and a free sample please call +44(0)1706 620600 or visit our website - www.tscswabs.co.uk

Issue 2: 24/05/2013



Technical Service
Consultants Ltd

Technical Service Consultants Ltd

Microbiology House, Fir Street, Heywood, Lancashire OL10 1NW, UK

T : +44 (0)1706 620600 F : +44 (0)1706 620445

E : sales@tscswabs.co.uk W : www.tscswabs.co.uk