

Press Release

Release date: December 2012

Quick and convenient solutions for environmental water sampling

Technical Service Consultants' offer a range of ready prepared sterile Environmental Sampling Containers with the addition of Sodium Thiosulphate for use in environmental and potable water sampling.


To make the process of neutralising chlorine in environmental and potable water samples simple, safe and cost effective, TSC have already included Sodium Thiosulphate to the sampling containers, so there is no additional preparation required by the end user.

500ml sampling containers are best used in routine microbiological testing of potable (drinking water) samples and 1000ml containers are ideal for sampling water from cooling towers. With a range of concentrations available, both containers are also suitable for sampling swimming pools, endoscopy rinse waters, hydrotherapy pools, renal dialysis water and dental line supplies.

Part of the reliable Enviroscreen Range, Environmental Sampling Containers are manufactured from durable polyethylene, have large labels to make identification of the sample site simple and straight forward and each sampling container has an irradiation indicator to confirm sterilisation.

All sampling containers have a 2 year shelf life and certificates of conformity can be downloaded directly from our website.

For further details and a free sample please visit www.tscswabs.co.uk or call +44(0)1706 620600.



Technical Service
Consultants Ltd

Technical Service Consultants Ltd

Microbiology House, Fir Street, Heywood, Lancashire OL10 1NW, UK
T : +44 (0)1706 620600 F : +44 (0)1706 620445
E : sales@tscswabs.co.uk W : www.tscswabs.co.uk