

# Press Release

## Technical Service Consultants offer simple and effective cleaning and sanitising validation systems.

Uncover hidden contamination by using The EnviroScreen range, offering customers a choice of ready to use sampling kits for routine microbiological and general hygiene monitoring of equipment and surfaces, ideal for validation of cleaning and sanitising programs.

### TSC Neutralising Buffer

In many cases a surface will contain traces of disinfectants and/or sanitisers. Left alone and over time, these cleaning fluid traces will not kill the microorganisms, but will just prevent any growth (bacteriostatic effect). TSC Neutralising Buffer will neutralise most common disinfectants, providing end-users with total testing assurance.

#### Quantitative Swab Sampling Kits

- High quality transport sampling kits
- Ready filled with TSC Neutralising Buffer
- Accurate volumes
- Choice of swab shafts and tips
- 2 year shelf life

#### PreMoist Hygiene Swab Kits

- Pre-moistened with TSC Neutralising Buffer, saving time & cost
- Entire sample can easily be analysed
- Choice of swab tips
- 12 month shelf life

#### Sampling Sponge Kits

- Pre-moistened with TSC Neutralising Buffer for convenience
- Fast & simple way to confirm the presence of living bacterial cells
- Cellulose sponge offers enhanced sample collection & release
- Sturdy breakpoint handle available for hard to reach areas
- 12 month shelf life
- Ideal monitoring solution

#### Dry Swab Sampling

- Wide selection of swab sizes, bud types and shafts
- With and without breakpoints for added convenience
- Ideal for general sampling applications of cleaned areas



All sampling devices conform to ISO 18593:2004 requirements and are gamma irradiated for assured sterility.  
For further details and samples please visit [www.tscswabs.co.uk](http://www.tscswabs.co.uk) or call +44(0)1706 620600

01/07/2014\_v2



**Technical Service**  
Consultants Ltd

#### Technical Service Consultants Ltd

Microbiology House, Fir Street, Heywood, Lancashire OL10 1NW, UK  
T : +44 (0)1706 620600 F : +44 (0)1706 620445  
E : [sales@tscswabs.co.uk](mailto:sales@tscswabs.co.uk) W : [www.tscswabs.co.uk](http://www.tscswabs.co.uk)