



Product Specification Sheet

SS-L03

SwabSURE Salmonellae for direct detection of Salmonellae by colour change

Product Code	SS-L03
Description	SwabSURE Salmonellae for direct detection of Salmonellae by colour change
Intended Use	Direct detection of Salmonellae within 18-24hours.
Components	75mm Polypropylene tube & label, polypropylene screw cap, Breakpoint viscose tip swab in 150mm polypropylene tube & Label.
Contents	De-ionised Water, Salmonellae Selective agents, Nutritive Broth, Dye, Tween 80, Dye, Lecithin, Sodium Thiosulphate, Sodium Thioglycollate, Buffers to maintain pH - 7.6±0.2.
Method of Sterilisation	Terminally Sterilised (media autoclaved, swab gamma irradiated)
Shelf Life	12 months from date of manufacture for media, 24 months for swabs.
Storage Conditions	4-8°C for Media and 4-25°C for swabs. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging
Precautions and Limitations	<p>1. Recommended use for the presumptive detection of Salmonellae from food contact & environmental surfaces.</p> <p>2. A positive colour change may arise due to the presence of excessively high numbers of competing organism species like <i>Citrobacter</i> spp. producing a false positive result. The occurrence of false positive results is an important indicator of suspect cleaning & sanitising practices due to the concentration of background flora & should flag a disinfectant review.</p> <p>3. Some uncommon strains (e.g. <i>S. arizonae</i>) may be H₂S negative and could be missed though are usually associated with animals (reptiles).</p>
Pack Size	100 vials and swabs in a box
Other	Quality control testing details on final product: i) Verification of sterility. ii) Target Organism recovery & Final colour

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 0 Date: 01.05.2022

ProductSpec_SS-L03