



Product Specification Sheet

BAG 3

100 Large Minigrip bags, folded & sealed inside 2 Large Minigrip bags - both outer bags tamper sealed with Irradiation indicators. 10 packs (1000) per box

Product Code	BAG 3
Description	100 Large Minigrip bags, folded & sealed inside 2 Large Minigrip bags - both outer bags tamper sealed with Irradiation indicators. 10 packs (1000) per box
Intended Use	Sample throughout. General laboratory use only
Components	Large Minigrips x110 (350mm x 230mm), Irradiation Sterility indicators x20. Recyclable Food grade LDPE.
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-53.4kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging
Precautions and Limitations	Not applicable.
Pack Size	1000 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 5 Date: 20.10.2020
ProductSpec_BAG3