



## Product Specification Sheet

### **TS/10-B**

Sterile swab with a wooden shaft and cottontip, 100 in sealed paper bag

Product Code	TS/10-B
Description	Sterile swab with a wooden shaft and cottontip, 100 in sealed paper bag
Intended Use	Sample throughout. General laboratory use only
Components	Cotton tip/woodshaft swab, paper bag, Non-toxic glue, Label & film.
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-53.4kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging
Precautions and Limitations	Not applicable.
Pack Size	8500 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 3 Date: 24.04.2020  
ProductSpec\_TS10-B