



## Product Specification Sheet

### **TS/14-C**

Sterile dry swab with viscose minitip and solid aluminium shaft, single in tube

Product Code	TS/14-C
Description	Sterile dry swab with viscose minitip and solid aluminium shaft, single in tube
Intended Use	Collection and transport of specimens
Components	Aluminium wire shaft/Cotton mini tip, orange standard cap, Polypropylene Tube, Non-toxic glue, Label & film.
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-53.4kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging
Precautions and Limitations	Not applicable.
Pack Size	100 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 3 Date: 24.04.2020  
ProductSpec\_TS14-C