



Product Specification Sheet

TS/15-2B

Poultry boot swab kit. 2 pairs premoistened in MRD all in resealable outer bag for sample transportation (includes overboots & gloves)

| | |
|-----------------------------|---|
| Product Code | TS/15-2B |
| Description | Poultry boot swab kit. 2 pairs premoistened in MRD all in resealable outer bag for sample transportation (includes overboots & gloves) |
| Intended Use | Poultry house litter testing for Salmonellae |
| Components | 4 boot swabs – 2 x Plastic overboots - 1 x Minigrip bag - 1 x gloves - 2 x Medical grade bag. |
| Contents | De-ionised Water – 1 Litre, Peptone – 1g, Sodium Chloride – 8.5g, pH – 7.0 ± 0.2. |
| Method of Sterilisation | Sterilised by gamma irradiation at dosage of 25-53.4kGy |
| Shelf Life | 12 months from date of manufacture |
| Storage Conditions | 4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging |
| Precautions and Limitations | 1. Peptone Saline Diluent (PSD) has been modified to enhance enrichment of target organisms. Further selective steps may be required. 2. Due to increased material absorbance over time, boots may vary in appearance as the liquid media settles into the deep |
| Pack Size | 25 in a box |
| Other | Quality control testing details on final product: i) Validation of irradiation sterility ii) Organism Recovery (Salmonella typhimurium). |

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 1 Date: 24.04.2020

ProductSpec_TS15-2B