



Product Specification Sheet

TS/15-BPW

Sterile blue sponge in easy open stomacher pouch, individually packed & dosed with Buffered Peptone Water

Product Code	TS/15-BPW
Description	Sterile blue sponge in easy open stomacher pouch, individually packed & dosed with Buffered Peptone Water
Intended Use	Sponge for pathogen monitoring in poultry houses
Components	Blue Cellulose Sponge, easy open stomacher pouch.
Contents	De-ionised Water, Bacteriological Peptone, Phosphate Buffers.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-53.4kGy
Shelf Life	12 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging
Precautions and Limitations	Buffered Peptone water (BPW) is a generalised media hence will recover all organisms. Further selective steps are required
Pack Size	200 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility. ii) Organism recovery (S.typhimurium)

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 4 Date: 24.04.2020

ProductSpec_TS15-BPW