



## Product Specification Sheet

### **TS/15-T40**

Sterile individually wrapped blue Plastic Template 10x10 cm with press 'hold' tab

Product Code	TS/15-T40
Description	Sterile individually wrapped blue Plastic Template 10x10 cm with press 'hold' tab
Intended Use	Defined area Sample Collection. General laboratory use only
Components	Blue Plastic, easy open stomacher pouch.
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-53.4kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging
Precautions and Limitations	Template can only be used once or cross over contamination will occur
Pack Size	40 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 5 Date: 24.04.2020

ProductSpec\_TS15-T40