



Product Specification Sheet

TS/15-U125

EnviroStik - Sponge with Straight Handle dosed in 10ml Neutralising Buffer in easy open stomacher pouch. No rivet design.

Product Code	TS/15-U125
Description	EnviroStik - Sponge with Straight Handle dosed in 10ml Neutralising Buffer in easy open stomacher pouch. No rivet design.
Intended Use	Environmental and hygiene monitoring
Components	Blue Cellulose Sponge, Blue plastic stick with built-in breakpoint, easy open stomacher pouch. Non toxic glue.
Contents	De-ionised Water, Tween 80, Lecithin, Sodium Thiosulphate, Sodium Thioglycollate, Buffers to maintain pH - 7.2 ± 0.2
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-53.4kGy
Shelf Life	12 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging
Precautions and Limitations	Neutralisation is limited to surface residues. Due to product characteristics a small precipitate may be observed which does not affect performance
Pack Size	125 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility ii) Product Specification iii) Organism recovery

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 0 Date: 12.09.2023
ProductSpec_TS15-U125