



Product Specification Sheet

TS/15-U:NaCl125

EnviroStik - Sponge with Straight Handle dosed in 10ml 0.9% Saline in easy open stomacher pouch. No rivet design.

Product Code	TS/15-U:NaCl125
Description	EnviroStik - Sponge with Straight Handle dosed in 10ml 0.9% Saline in easy open stomacher pouch. No rivet design.
Intended Use	Sampling of meat carcasses for microbiological analysis
Components	Blue Cellulose Sponge, Blue plastic stick with built-in breakpoint, easy open stomacher pouch. Non toxic glue.
Contents	De-ionised Water – 1 Litre, Sodium Chloride –9g.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-53.4kGy
Shelf Life	12 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging
Precautions and Limitations	Not applicable.
Pack Size	125 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility ii) Product Specification iii) Organism recovery

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 0 Date: 12.09.2023
ProductSpec_TS15-U:NaCl125