



## Product Specification Sheet

### **TS/16-A250**

Sterile dry treated paper shaft, Cotton mini tip Swab, single in tube.

Product Code	TS/16-A250
Description	Sterile dry treated paper shaft, Cotton mini tip Swab, single in tube.
Intended Use	Culture swab; recommended for urethral sampling
Components	plasticised paper shaft/cotton mini tip swab, orange cap, Polypropylene Tube, Non-toxic glue, Label.
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-53.4kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging
Precautions and Limitations	Not applicable.
Pack Size	250 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 4 Date: 24.04.2020  
ProductSpec\_TS16-A250