



## Product Specification Sheet

### **TS/16-S**

Sterile dry swab with cotton mini tip and treated paper shaft, breakpoint at 45mm, single in peel pouch.

Product Code	TS/16-S
Description	Sterile dry swab with cotton mini tip and treated paper shaft, breakpoint at 45mm, single in peel pouch.
Intended Use	Culture swab; recommended for urethral sampling
Components	Plasticised paper shaft/cotton mini tip (breakpoint at 45mm), peel pouch film.
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-53.4kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging
Precautions and Limitations	Recovery after 24 hours will be unpredictable and target organisms will be lost.
Pack Size	1000 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 5 Date: 24.04.2020  
ProductSpec\_TS16-S