



Product Specification Sheet

TS/165

Sterile 4000 ml Non-sterile Empty Sampling Container for 24 hour Urine Analysis

Product Code	TS/165
Description	Sterile 4000 ml Non-sterile Empty Sampling Container for 24 hour Urine Analysis
Intended Use	Multipurpose Use
Components	HDPE Bottle, HDPE Cap, Label.
Contents	Not applicable as supplied empty.
Method of Sterilisation	N/A - supplied non-sterile
Shelf Life	N/A
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging
Precautions and Limitations	The bottle is made of High Density Polyethylene and is only suitable for water and certain solvent/chemical treatments. Please check compatibility before use.
Pack Size	50 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 3 Date: 14.11.2022
ProductSpec_TS165