



Product Specification Sheet

TS/17-H250

Swab with a flexible polypropylene shaft tipped with pure viscose, in polypropylene outer tube.

Product Code	TS/17-H250
Description	Swab with a flexible polypropylene shaft tipped with pure viscose, in polypropylene outer tube.
Intended Use	Culture swab; Hygiene sampling & general Laboratory use
Components	Polypropylene Tube, HDPE cap, polypropylene/viscose tip swab, Non-toxic glue, Label.
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-53.4kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging
Precautions and Limitations	Not applicable.
Pack Size	250 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 1 Date: 24.04.2020

ProductSpec_TS17-H250