



Product Specification Sheet

TS/18

3 inch cotton tipped swab with a flexible polypropylene shaft, 100 swabs in Minigrip bag.

Product Code	TS/18
Description	3 inch cotton tipped swab with a flexible polypropylene shaft, 100 swabs in Minigrip bag.
Intended Use	Ideal for general swabbing where a small shaft is required.
Components	Polypropylene/Cotton tip swab, Minigrip bag.
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-53.4kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging
Precautions and Limitations	Not applicable.
Pack Size	5000 in a box (50 packs of 100)
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 3 Date: 24.04.2020

ProductSpec_TS18