



## Product Specification Sheet

### **TS/19-H250**

Blue polystyrene Breakpoint shaft @45mm, Viscose tip Swab, in Tube

Product Code	TS/19-H250
Description	Blue polystyrene Breakpoint shaft @45mm, Viscose tip Swab, in Tube
Intended Use	Culture swab; general Laboratory use - Food & Hygiene use
Components	Polypropylene Tube, Blue Standard Cap, Polystyrene/Viscose swab, Non-toxic glue, Label.
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-53.4kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging
Precautions and Limitations	Not applicable.
Pack Size	250 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 3 Date: 24.04.2020  
ProductSpec\_TS19-H250