



Product Specification Sheet

TS/19-S

Polystyrene shaft with breakpoint @45mm, Viscose tip swab in peel pouch.

Product Code	TS/19-S
Description	Polystyrene shaft with breakpoint @45mm, Viscose tip swab in peel pouch.
Intended Use	Culture swab; general Laboratory use
Components	Polystyrene/Viscose tip swab (BP@45mm), peel pouch film.
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-53.4kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging
Precautions and Limitations	Not applicable.
Pack Size	1000 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 4 Date: 24.04.2020
ProductSpec_TS19-S