



Product Specification Sheet

TS/25

Plasticised Paper Shaft/Cotton MiniTip (Male), Polystyrene Shaft/Viscose Tip(Female), in Peel Pouch

Product Code	TS/25
Description	Plasticised Paper Shaft/Cotton MiniTip (Male), Polystyrene Shaft/Viscose Tip(Female), in Peel Pouch
Intended Use	Collection and transport of specimens including Chlamydia spp. - dual swab pack
Components	Plasticised paper/Cotton mini tip swab (Male), Polystyrene/Viscose Swab (Female), Peel Pouch film.
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-53.4kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging
Precautions and Limitations	Not applicable.
Pack Size	1000 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 5 Date: 24.04.2020
ProductSpec_TS25