



Product Specification Sheet

TS/30-C

Sterile L shaped Spreaders (firm), 25 in a sealed minigrip bag

Product Code	TS/30-C
Description	Sterile L shaped Spreaders (firm), 25 in a sealed minigrip bag
Intended Use	Sample throughout or Protect. General laboratory use only
Components	L shaped Plastic (firm), Minigrip bag.
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-53.4kGy
Shelf Life	36 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging
Precautions and Limitations	Sterility is only guaranteed whilst bag is sealed. Once opened, any unused items should be discarded
Pack Size	1000 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 5 Date: 24.04.2020
ProductSpec_TS30-C