



## Product Specification Sheet

### **TS/4-D**

Sterile 1µl Yellow Culture Loops (soft), 20 in the sealed minigrip bag

Product Code	TS/4-D
Description	Sterile 1µl Yellow Culture Loops (soft), 20 in the sealed minigrip bag
Intended Use	Sample throughout. General laboratory use only
Components	1 µl Plastic Yellow Loops (soft), Minigrip bag.
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 20-53.4kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging
Precautions and Limitations	Laboratory use only
Pack Size	1000 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 5 Date: 24.04.2020  
ProductSpec\_TS4-D