



## Product Specification Sheet

### **TS/5-18A**

Probact Amies transport medium with charcoal in 75mm tube with black cap  
75mm polystyrene shaft/pure cotton swab in peel pouch

Product Code	TS/5-18A
Description	Probact Amies transport medium with charcoal in 75mm tube with black cap 75mm polystyrene shaft/pure cotton swab in peel pouch
Intended Use	Bacteriological samples only. Collection and transport of clinical specimens
Components	Polypropylene Tube (80 mm), Blue HDPE Standard Cap, Black HDPE Standard Cap, Polystyrene/Viscose swab (shaft length 80mm), Non-toxic glue, Label, peel pouch film.
Contents	De-ionised water, Amies Charcoal Transport Medium, pH – $7.15 \pm 0.2$ .
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-53.4kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging
Precautions and Limitations	Recovery after 48 hours will be unpredictable and target organisms will be lost. Fatty acid sensitive strains will give decreased recovery or be stressed in the presence of cotton
Pack Size	100 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility ii) Organism recovery (N.gonorrhoeae, H.influenzae)

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.