



Product Specification Sheet

TS/5-9

Probact Amies transport medium without charcoal in tube with blue cap polystyrene/pure viscose 45mm breakpoint swab in peel pouch

Product Code	TS/5-9
Description	Probact Amies transport medium without charcoal in tube with blue cap polystyrene/pure viscose 45mm breakpoint swab in peel pouch
Intended Use	Bacteriological samples only. Collection and transport of clinical specimens
Components	Polypropylene Tube, Blue HDPE Standard Cap x2, Polystyrene/Viscose swab, Non-toxic glue, Label & peel pouch film.
Contents	De-ionised water, Amies without Charcoal Transport Medium, pH - 7.15±0.2.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-53.4kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging
Precautions and Limitations	Recovery after 24 hours will be unpredictable and target organisms will be lost. Fatty acid sensitive strains will give decreased recovery or be stressed in the presence of cotton
Pack Size	100 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility ii) Organism recovery

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 1 Date: 24.04.2020
ProductSpec_TS5-9