



## Product Specification Sheet

### **TS/6-S6:250**

Sterile woodshaft/sharpened cotton swab sealed in outer tube

|                             |  |
|-----------------------------|--|
| Product Code                | TS/6-S6:250  |
| Description                 | Sterile woodshaft/sharpened cotton swab sealed in outer tube   |
| Intended Use                | Culture swab; general laboratory use; forensic or precision work   |
| Components                  | Polypropylene Tube, Blue HDPE Standard Cap, Woodshaft/Tapered Cotton tip, Non-toxic glue, Label.                     |
| Contents                    | Not applicable.  |
| Method of Sterilisation     | Sterilised by gamma irradiation at dosage of 25-53.4kGy  |
| Shelf Life                  | 24 months from date of manufacture   |
| Storage Conditions          | 4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging |
| Precautions and Limitations | Not applicable.  |
| Pack Size                   | 250 in a box   |
| Other                       | Quality control testing details on final product: i) Validation of irradiation sterility.                            |

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 3 Date: 24.05.2020  
ProductSpec\_TS6-S6:250