



Product Specification Sheet

TS/9-B

Sterile Woodshaft/cotton swab, 5 in peel pouch

Product Code	TS/9-B
Description	Sterile Woodshaft/cotton swab, 5 in peel pouch
Intended Use	Culture swab & general laboratory use
Components	Woodshaft/Cotton swab, Peel pouch film.
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-53.4kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging
Precautions and Limitations	Once opened, pack must be used aseptically & discard any swabs left over
Pack Size	1000 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 5 Date: 24.04.2020
ProductSpec_TS9-B